

Delivered within 5 working days*



in partnership with



Credalast Circular Knit Made to Measure Compression Stockings. Available on Drug Tariff

Order form for Credalast Velvet and Credalast Classic Made to Measure Stockings

It is essential that ALL Parts 1, 2, 3 & 4 are FULLY COMPLETED, otherwise your order cannot be processed

PART 1 (BLOCK CAPITALS PLEASE)

Name of pharmacy: _____

Account/Branch number: (if known) _____

Pharmacy address: _____

Postcode: _____

Pharmacy telephone number: _____


E-mail address: _____

Pharmacy contact: _____

Patient's name: _____

Prescription number: _____

Patient's previous reference no.: (if known) SM/ _____



PART 2 Please complete ALL sections and tick the appropriate box:

Section One: Choose which range of stocking
 SEE OVERLEAF FOR CHOICE **Credalast Velvet MTM**
 Credalast Classic MTM

Section Two: Compression: Class I (14-17mm Hg)
 Class II (18-24mm Hg)
 Class III (25-35mm Hg)

Section Three: Shade: Natural
 Black

Section Four: Foot: Closed Toe (CA/LZ)
 Open Toe (CA/LA)

Section Five: Style: Thigh Length
 Start measurements at CG/LG
 Below Knee
 Start measurements at CD/LD

PART 3 QUANTITY

| | Left Leg | Right Leg |
|--|----------------------|----------------------|
| (MUST BE COMPLETED WITH NUMBERS OR ZERO IN EACH BOX) | <input type="text"/> | <input type="text"/> |

NOTES

PART 4 CA measurement must be completed for closed and open toe

| C = CIRCUMFERENCE | Left cm | Right cm |
|---|---------|----------|
| CG (Thigh) | | |
| CF (Mid Thigh) | | |
| CD (Below Knee) | | |
| CC (Calf) | | |
| CB (Ankle) | | |
| CY (Circumference around the heel) | | |
| CA (Circumference of foot measured at base of toes) | | |
| L = LENGTH | Left cm | Right cm |
| LG | | |
| LF | | |
| LD | | |
| LC | | |
| LB | | |
| LA (Open toe) | | |
| LZ (Closed toe) | | |

*Items will be dispatched within 5 working days if all details are correct and subject to factors beyond our control.

ALL COMPLETED ORDERS INCLUDING REPEAT REQUESTS TO BE E-MAILED TO MTM@CREDENHILL.CO.UK OR FAXED ON **0115 944 0437**
 WE REGRET WE ARE UNABLE TO PROCESS TELEPHONE ORDERS

To download additional forms go to www.credalast.co.uk
 Form designed to be faxed to Credenhill by the organisation dispensing the stockings.

Submission of this form confirms that the patient's consent has been obtained to submit these details to the manufacturer in order to produce the patient's made to measure stockings.

All measurements should be taken on the bare leg whilst standing or following a period of rest. The manufacturer cannot be held responsible for incorrect measurements.

Delivered within 5 working days*

Delivered within 5 working days*




in partnership with



Credalast Circular Knit Made to Measure Compression Stockings. Available on Drug Tariff

Guidance notes for Credalast Velvet and Credalast Classic Made to Measure Stockings



Credalast stockings are available on **drug tariff FP10**

Please ensure **Parts 1-4** are completed overleaf along with the **prescription** (each option below must be listed on the prescription)

- Velvet or Classic • Compression Class • Colour • Open or Closed Toe
- Below Knee or Thigh Length • Quantity for each leg

If you are a **Pharmacist**: email measurement form only to mtm@credenhill.co.uk (scans or smart phone pictures excepted)

If you are a **Patient** or **Nurse Prescriber**: Send MTM and prescription to FREEPOST RTHL-YUAG-GCST Daylong, Credenhill Ltd, Unit 10, Cossall Industrial Estate ILKESTON, DE7 5UG

Alternatively fax the forms on **0115 944 0437**. Hosiery will normally arrive **within 5 working days** of Credenhill receiving the form. If you need any further help please ring us on **0115 932 0144**, our experienced team will be able to help.


Credenhill currently offers a choice of two stockings, **Credalast® Velvet** and **Credalast® Classic** both developed from premium Swiss stocking manufacturing technology.

Credalast® Velvet

The first prescription Made to Measure stocking available with the Tactel® Climate Effect which wicks away moisture and improves air permeability, helping to keep an optimum skin balance regardless of activity.

Credalast® Classic

This traditional stocking offers a more conventional stiffness and style which suits many patients and is still available.



Application Aids

Thigh Some users may find it difficult to apply compression hosiery. For this reason we stock a range of highly useful application aids such as Easy Slide designed to improve your wearer experience.

Accessories

Thigh length stockings will need suspending as they are not manufactured with a self supporting top band.

We offer Credalast Suspender Belt available on drug tariff FP10 or Credalast Y Shaped Suspenders. We can also provide Coltex – a skin adhesive designed to hold compression hosiery in place.

Credalast Classic® and Credalast Velvet® are registered trade marks. Tactel® is a registered trade mark of Invista.

© 2017. Credenhill. All Rights Reserved

CR001 V.8-09/17

Credenhill Limited, 10 Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UG.

Telephone: 0115 932 0144

Email: mtm@credenhill.co.uk